**GR
EN**

IFU rev2..0_07/2018

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**MOBIAK**[®]
MEDICAL • GASES • FIREFIGHTING est. 1977

INTRODUCTION

Thank you for purchasing a wheelchair by MOBIAK A.E product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance, In case of doubts, we recommend to contact the dealer who can help and advice you properly.

Note: check all parts for shipping damage and test before using. In case of damage, do not use the device. Contact your dealer further instruction.

PRODUCT DESCRIPTION

These wheelchairs have been designed for physically or mentally handicapped people with limited or no personal mobility only for transport.

This user's manual is compiled from the latest specifications and product information available at the time of publication. We reserve the right to make changes as they become necessary. Any changes to our products may cause slight variations between the illustrations and explanations in this manual and the product you have purchased.

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Meets the requirements of the Medical Device Directive: 93/42/EEC

We declare the conformity of the above product to the provisions of the Medical Equipment Directive 93/42/EEC and the Greek National Legislation ΔΥ8δ/Γ.Π.οικ.130648. ΦΕΚ 2198Β/02-10-2009

We follow the ISO 13485, Quality Management System for Medical Devices & ISO 9001 Quality Management System

OTHER INSTRUCTIONS

- The wheelchairs have a brake on the back wheels, which should be secured when the wheelchair is not moving.
- DO NOT push on the footrests to sit or stand up, DANGER OF OVERTHROW.
- The footrests are used by the user when they are seated.
- When standing up, you must rely on both handles equally. Supporting only one hand (and unattended) can cause injury.
- DO NOT leave the wheelchair wet.
- DO NOT leave children near the wheelchair unattended.
- DO NOT touch moving parts while moving DANGER OF INJURY.

For individual pieces only refer to the MOBIK A.E list

Getting in the wheelchair

- Position the wheelchair as close as possible to you
- Lean on the frame with arms and slowly go down towards the seat
- Engage wheel locks

Getting out of the wheelchair

- Position the wheelchair as close as possible along side the seat to which you are transferring
- Shift body weight into seat with transfer

MAINTENANCE

Periodic review of the procedures below will help to maintain the life and effectiveness of your wheelchair.

WEEKLY

- Check the wheel pressure by following the instructions on the wheels.
- Check that the central bolt of the transverse arm is tight and in good condition.
- Clean the wallpaper with soapy water or sufficient products.
- Clean the metal parts of the chair with a dry cloth or with special products for chrome surfaces or painted surfaces.
- Check the wheel bolts and adjust if necessary.
- Check for damage to front and rear wheels.

QUARTERLY

- Check and adjust if necessary: wheel forks, pulling wheels, propellers.
- Grease or oil all components that are subject to friction or impending wear.
- Tighten the screws and nuts.
- General revision of the chair's status and repair to whatever is necessary.

ANNUALLY

We recommend a full review once a year to a qualified service provider.

While moving with your wheelchair, do not forget !!

Always use lifts and ramps. If these are not available, you can reach your destination with the help of two escorts. They must catch the wheelchair in securely seated areas. A wheelchair should never be lifted by a passenger. It should be pushed only.

CLEANING

In addition to maintenance and inspection, both the frame and the wallpaper can be thoroughly cleaned.

- To maintain its appearance, the wallpaper must be cleaned regularly.
- Dirt can be removed with a sponge and a mild detergent with warm water, sweeping in the direction of texture.
- Rinse with clean water before letting it dry.
- Do not place the wallpaper in a tumble dryer. If the dirt is deep, the surface should be cleaned with a soft brush.
- The frame should be cleaned regularly with a mild detergent.
- The wheelchairs with leather seats must be cleaned with dry cloth.
- Clean the wheels with warm water and mild soap, Pay attention in drying the wheels , it's dangerous use the wheelchair with wet wheels.

SAFETY ADVICE

What to check before use

- Quick-release axles (if any) on the rear wheels.
- Tires, tire pressures and wheel locks. The wheel locks will not work properly if the pressure is too low or the gap between the tire and the wheel locks is too long.
- The maximum user weight is 120 kg. This wheelchair should only be used to carry one person each time on the seat. Any other use beyond this limit violates the security regulations.
- Users should not sit in the wheelchair while driving on any means of transport.
- When vehicles are in motion, unoccupied wheelchairs should be secured by appropriate means. It is likely that traps or uneven terrain could lead the wheelchair to skip, especially when driving uphill or downhill.
- When driving in a step or tilting forward, the body must lean forward.

IMPORTANT!

- The construction of this wheelchair are designed to provide maximum safety.
- The international safety standards currently in force have either been fulfilled or exceeded. However, users may be at risk by using their wheelchairs unsuitable.
- For your own safety, the following rules must be observed.
- Nonspecific or incorrect adjustment work can increase the risk of accidents.
- Always wear light-colored clothes when moving in dark places so that you can be reached more easily from others.
- Be careful during your first course in this wheelchair.
- Meet your wheelchair.
- Make use of your wheelchair on level ground as well as uphill and downhill.
- Do not attempt to reach objects if you need to lift them from the floor reaching between your knees.
- Do not attempt to reach objects if you need to move on to the seat.
- Do not use on roads or motorways.
- Do not attempt to climb or descend with water, ice or oil membrane.
- Do not use unauthorized parts, accessories or adapters other than those approved by Mobiak.
- Do not stand in the wheelchair.
- Do not use the sideboard as a platform when entering or exiting the wheelchair.

WARNING

In case of a problem , do not use the wheelchair until the problem will be solved.

- Do not use the product for a purpose not indicated in this manual.
- MOBIAK declines all responsibilities for any consequences resulting from an incorrect use of this product.
- The manufacturer reserves the right to change the information contained inn this document without notice.

- Advise with physician or therapist to determinate the correct adjustment and the correct use of the wheelchair.
- Do not allow children to play on the wheelchair.
- Do not attempt to lift the wheelchair by any removable parts.
- Keep the packed wheelchair away from heat sources.
- Getting in and out of the wheelchair may be performed independently , provided you have an adequate mobility and upper body strength.

Protecting the environment

Should you find one day that your appliance needs to be replaced or if it is of no further use to you, think of the protection of the environment:

1) Do not dispose your appliance along with the rest of the public waste (this is also the meaning of the shown recycling sign).

2) Contact your Public Authorities and they will instruct you of the Recycling centers to which your appliances must be disposed.

3) Correct disposal of your appliance helps the protection of the environment as well as the recycling of the appliance components.

4) The electrical appliances due to their construction materials, if not correctly disposed, may lead to environmental and furthermore health hazards.

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EN. First of all we check the contents of the package. Place the wheelchair as shown in the picture.

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EN. By pushing the bars on the side of the seat with one movement, the wheelchair opens and comes into use.

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EN. Dragging the handles where they are at the same point, with a simple move, the wheelchair comes into a state of safekeeping.

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EN. In case your wheelchair is accompanied by a pillow (**0808480 - 0808481 - 0808482**), place the Velcro side on the strips on the seat.

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EN. Move the armrest forward clockwise until the bottom of the seat meets the front of the seat.

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EN. At the point where the armrest meets the seat, there is a reset lever. Pull it back and push the armrest down to lock it.

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EN. Release the lever to return to snap the hole. Now the armrest is locked and secure. The reverse process is released.

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EN. Hold the footrest with your hand and place it in the way shown. The pin has to pass through its hole.

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EN. Holding the footrest firmly, place it so that the parallel with the hole part touches and touches the vertical tube where it ends up on the front wheel.

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EN. Once the curved section of the bottom is tangential to the vertical tube, rotate the footrest counterclockwise until the lock switch is closed and locked.

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EN. Use the straps included to make sure the patient's feet are properly positioned.

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EN. These straps are fitted with a Velcro fastener and placed vertically on the feet.

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EN. Our wheelchair must be in the form shown, so it is ready for use. Below we will see some details where we need to be noticed.

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EN. By pushing forward the lever to the right and left of the AA, in front of the rear wheel, the brake of the respective wheel is activated. If the lever moves sufficiently forward, the brake locks for prolonged stability of the AA. With reverse lever movement, the brake is released and the wheelchair can move again.

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EN. The user must move the AA through the drive shafts where they are mounted on the wheel from the outside of the wheel. Otherwise, there is a risk that the user may be injured when moving the wheelchair.

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EN. In the case where the AA includes an overturning system, remove the cap from the pipes at the bottom of the wheelchair, parallel to the rear wheels.

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EN. Insert rolling overturning until the hole in the trunk intersects the hole in the tube where the cap was removed.

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

EN. Install the lock pin where it is included so that it passes through the intersecting holes and appears on the other side of the pipe. By inverting the process we can remove the rollover avoidance system.

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EN. Codes with commode, check if all spare parts are included in packaging and place it in the corresponding sockets.

WARRANTY:

Your wheelchair is guaranteed for **2 years**. No service or spare part is covered unless it has been inspected by an authorized service or the importer / manufacturer of the product. The warranty covers factory failures & NO damage from misuse - maintenance of the carriage. Shipping costs for goods of warranty are to be paid by the customer or by dealer.

Warranty Claim form

Buyers Data

Name :

Date of purchase:

Phone:

Invoice number:.....

Resellers Data

Name :

Address :.....

Phone:

Reseller Stamp

Warranty Exceptions

- Belts, straps, pillows and other subject to daily contact with the user
- Floor Standings ..
- Wheel Tires.
- Any damage has been caused by misuse, dropping or hitting the wheelchair
- Any damage caused by poor maintenance, transportation or storage of the wheelchair.
- Any normal wear and tear caused by the use of the wheelchair .