

Mask Quick Reference Guide



Mask Quick Reference Guide



Nasal Masks and Spare Parts

F&P Eson™

F&P Zest™ Q

F&P Lady Zest™ Q

Zest™

FlexiFit™ 407

FlexiFit™ 406 Petite

FlexiFit™ 405

Full Face Masks and Spare Parts

F&P Simplus™

Forma™

FlexiFit™ 431

FlexiFit™ 432

Nasal Pillows Masks and Spare Parts

F&P Pilairo™ Q

Opus™ 360

Oral Mask and Spare Parts

Oracle™

F&P Eson™



The F&P Eson is designed to perform in tune with you and your needs. The three key components, the *RollFit™ Seal*, *ErgoFit Headgear* and *Easy Frame*, work in harmony. These components have been designed to deliver the comfort, seal and easy use that Fisher & Paykel Healthcare masks are known for.

- **RollFit Seal**

As the name suggests, the one-piece seal ‘rolls’ back and forth on the bridge of the nose to adjust automatically. In doing so, the RollFit Technology reduces pressure on the bridge of the nose.

- **ErgoFit Headgear**






This breathable, ergonomically designed headgear has been designed to self-locate high on the rear of the head, allowing for head movement (sideways, up and down).

- **Easy Frame**

This low-profile frame is stable, durable and small, and ensures a clear line of sight. The one frame fits all three seal sizes and has an Easy-Clip frame attachment to aid assembly after cleaning.

Product Code	Description	Quantity	HCPCS Code
400449 F&P Eson Mask Small	F&P Eson Nasal Mask (Complete System) includes:		A7034 + A7035
400450 F&P Eson Mask Medium			
400451 F&P Eson Mask Large			
	• RollFit Seal	1	
	• ErgoFit Headgear	1	
	• Diffuser Material plus spares	6	
	• Q Cover	1	

The F&P Eson Nasal Mask spare parts are listed below.

Product Code	Description	HCPCS Code	
400HC567* 400HC568*	Headgear – Small Headgear – Medium/Large	A7035	
400HC126* 400HC127* 400HC128*	Seal – Small Seal – Medium Seal – Large	A7032	
400HC228	Diffuser – 10 Pack + 1 Q Cover		
400HC569	Buckle & Clips		
400HC570	Swivel		
400HC571* 400HC572* 400HC573*	Mask Kit – Small (no headgear) Mask Kit – Medium (no headgear) Mask Kit – Large (no headgear)	A7034	

* Reimbursable replacement parts

F&P Zest™ Q

Product Code	Description	HCPCS Code	
400444 400445 400446	F&P Zest Q Petite F&P Zest Q F&P Zest Q Plus Complete system includes: <ul style="list-style-type: none"> • Easy-Clip Silicone Seal • FlexiFoam™ Cushion • Stretchgear™ Headgear • Glider™ Strap • Exhaust Flow Diffuser Cover • Exhaust Flow Diffuser Material – 5 pack + 1 	A7034 + A7035	
400HC557* 400HC542* 400HC558*	FlexiFoam Cushion and Silicone Seal – Petite FlexiFoam Cushion and Silicone Seal FlexiFoam Cushion and Silicone Seal – Plus	A7032	
400HC314* 400HC316*	Headgear Headgear – Petite	A7035	
400HC225	Exhaust Flow Diffuser Cover		
400HC226	Exhaust Flow Diffuser Material – 10 pack		
400HC227	Elbow (includes Swivel, Cover and Diffuser Material)		
900HC418	Glider Strap		
400HC576* 400HC577* 400HC578*	Mask Kit – Petite (no headgear) Mask Kit (no headgear) Mask Kit – Plus (no headgear)	A7034	

* Reimbursable replacement parts

F&P Lady Zest™ Q

Product Code	Description	HCPCS Code	
400447 400448	F&P Lady Zest Q Petite F&P Lady Zest Q Complete system includes: <ul style="list-style-type: none"> • Easy-Clip Silicone Seal • FlexiFoam™ Cushion • Rose Headgear • Glider™ Strap • Exhaust Flow Diffuser Cover • Exhaust Flow Diffuser Material – 5 pack + 1 	A7034 + A7035	
400HC542* 400HC557*	FlexiFoam Cushion and Silicone Seal FlexiFoam Cushion and Silicone Seal – Petite	A7032	
400HC325*	Headgear – Rose	A7035	
400HC225	Exhaust Flow Diffuser Cover		
400HC226	Exhaust Flow Diffuser Material – 10 pack		
400HC227	Elbow (includes Swivel, Cover and Diffuser Material)		
900HC418	Glider Strap		
400HC576* 400HC577*	Mask Kit – Petite (no headgear) Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

Zest™

Product Code	Description	HCPCS Code	
400439 400440 400441	Zest Petite Zest Zest Plus Complete system includes: <ul style="list-style-type: none"> • Easy-Clip Silicone Seal • FlexiFoam Cushion • Stretchgear Headgear • Glider Strap 	A7034 + A7035	
400HC542* 400HC557* 400HC558*	FlexiFoam Cushion and Silicone Seal FlexiFoam Cushion and Silicone Seal – Petite FlexiFoam Cushion and Silicone Seal – Plus	A7032	
400HC314* 400HC316*	Headgear Headgear – Petite	A7035	
400HC221	Elbow plus Swivel		
900HC418	Glider Strap		
400HC543* 400HC559* 400HC560*	Mask Kit (no headgear) Mask Kit – Petite (no headgear) Mask Kit – Plus (no headgear)	A7034	

* Reimbursable replacement parts

FlexiFit™ 407

Product Code	Description	HCPCS Code	
HC407	FlexiFit 407 Nasal Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (one size) • Foam Cushion • Headgear • Glider Strap 	A7034 + A7035	
400HC501*	Foam Cushion and Silicone Seal	A7032	
400HC301*	Headgear (one size)	A7035	
400HC202	Exhaust Flow Elbow		
900HC418	Glider Strap		
400HC502*	Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

FlexiFit™ 406

Product Code	Description	HCPCS Code	
HC406	FlexiFit 406 Petite Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (one size) • Foam Cushion • Headgear • Glider Strap 	A7034 + A7035	
400HC509*	Foam Cushion and Silicone Seal	A7032	
400HC307*	Headgear (Petite)	A7035	
400HC202	Exhaust Flow Elbow		
900HC418	Glider Strap		
400HC510*	Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

FlexiFit™ 405

Product Code	Description	HCPCS Code	
HC405	FlexiFit 405 Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small) • Silicone Seal (Large) • Foam Cushion • Headgear (one size fits all) • Exhaust Flow Diffuser Cap plus Spare • Exhaust Flow Diffuser Material plus Spares • Glider Strap 	A7034 + A7035	
900HC427* 900HC428*	Foam Cushion and Silicone Seal – Small Foam Cushion and Silicone Seal – Large	A7032	
900HC429*	Headgear (one size)	A7035	
900HC439	Exhaust Flow Diffuser Material – 10 Pack		
900HC440	Exhaust Flow Diffuser Cap plus Spare		
900HC441	Exhaust Flow Elbow and Diffuser Material with Spares		
900HC418	Glider Strap		
900HC406*	Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

F&P Simplus™



When engineering the F&P Simplus our designers set out to create a mask that revolutionized Full Face comfort, seal and easy use. The F&P Simplus is the result of intensive research and development, many prototypes and thousands of design hours. We consulted patients, equipment providers and clinicians from around the world to develop a mask that simply fits + performs.

- **RollFit Seal**

As the name suggests, the one-piece seal “rolls” back and forth on the bridge of the nose to adjust automatically. In doing so, the **RollFit Technology** reduces pressure on the bridge of the nose.

- **ErgoForm Headgear**






This breathable, ergonomically designed headgear self-locates high on the rear of the head, allowing for unrestricted head movement (sideways, up and down).

- **Easy Frame**

This low-profile frame is stable, durable and small, and ensures a clear line of sight. The one frame fits all three seal sizes and has an **Easy-Clip Frame Attachment** to aid assembly after cleaning.






Product Code	Description	Quantity	HCPCS Code
400475 F&P Simplus Small	Simplus Full Face Mask (Complete System) includes:		A7030 + A7035
400476 F&P Simplus Medium			
400477 F&P Simplus Large			
	• RollFit™ Seal	1	
	• ErgoForm™ Headgear	1	
	• Easy Frame	1	
	• Swivel	1	

The F&P Simplus Full Face Mask main spare parts are listed below.

Product Code	Description	HCPCS Code	
400HC579* 400HC580* 400HC581*	F&P Simplus RollFit™ Seal - Small F&P Simplus RollFit™ Seal - Medium F&P Simplus RollFit™ Seal - Large	A7031	
400HC582* 400HC583*	ErgoForm Headgear - Small ErgoForm Headgear – Medium/Large	A7035	
400HC587	F&P Simplus Headgear Clips + Buckle		
400HC570	F&P Gray Swivel		
400HC584* 400HC585* 400HC586*	F&P Simplus Mask – No Headgear Small F&P Simplus Mask – No Headgear Medium F&P Simplus Mask – No Headgear Large	A7030	

* Reimbursable replacement parts

Forma™

Product Code	Description	HCPCS Code	
400470 400471 400473	Forma – Small Forma – Medium and Large Forma – X-Large Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small or Medium and Large or X-Large) • FlexiFoam Cushion (Small or Medium and Large or X-Large) • Headgear with Easy-Release Cord • Glider Strap 	A7030 + A7035	
400HC121* 400HC122* 400HC123* 400HC124*	FlexiFoam Cushion and Silicone Seal – Small FlexiFoam Cushion and Silicone Seal – Medium FlexiFoam Cushion and Silicone Seal – Large FlexiFoam Cushion and Silicone Seal – X-Large	A7031	
400HC315*	Headgear with Crown Strap (one size)	A7035	
400HC203	Elbow and Vented Non-Rebreathing Valve		
400HC552* 400HC553* 400HC554* 400HC555*	Small Mask Kit (no headgear) Medium Mask Kit (no headgear) Large Mask Kit (no headgear) X-Large Mask Kit (no headgear)	A7030	

* Reimbursable replacement parts

FlexiFit™ 431

Product Code	Description	HCPCS Code	
HC431	FlexiFit 431 Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small) • Silicone Seal (Medium) • Silicone Seal (Large) • Headgear with Easy-Release Cord • Glider Strap 	A7030 + A7035	
400HC112* 400HC104* 400HC103*	Silicone Seal – Small Silicone Seal – Medium Silicone Seal – Large	A7031	
400HC302*	Headgear with Crown Strap (one size)	A7035	
400HC203	Elbow and Vented Non-Rebreathing Valve		
400HC503*	Mask Kit (no headgear)	A7030	

* Reimbursable replacement parts

FlexiFit™ 432

Product Code	Description	HCPCS Code	
HC432S HC432M HC432L HC432XL	FlexiFit 432 – Small FlexiFit 432 – Medium FlexiFit 432 – Large FlexiFit 432 – X-Large Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small or Medium or Large or X-Large) • Foam Cushion (Small or Medium or Large or X-Large) • Headgear with Easy-Release Cord • Glider Strap 	A7030 + A7035	
400HC113* 400HC114* 400HC115* 400HC119*	Foam Cushion and Silicone Seal – Small Foam Cushion and Silicone Seal – Medium Foam Cushion and Silicone Seal – Large Foam Cushion and Silicone Seal – X-Large	A7031	
400HC309*	Headgear with Crown Strap (one size)	A7035	
400HC203	Elbow and Vented Non-Rebreathing Valve		
400HC519* 400HC521* 400HC522* 400HC541*	Small Mask Kit (no headgear) Medium Mask Kit (no headgear) Large Mask Kit (no headgear) X-Large Mask Kit (no headgear)	A7030	

* Reimbursable replacement parts

F&P Pilairo™ Q



Weighing in as one of the lightest nasal pillows masks (1.94 ounces/55 grams¹), the F&P Pilairo Q has been further developed and re-engineered to improve patient satisfaction. The F&P Pilairo Q has just three key components and two minimalist headgear options.

- **AirPillow™ Seal**

The AirPillow Seal is 3D contoured with fine Variable Thickness Silicone (VTS™) and uses CPAP pressure to self-inflate to create an effective yet gentle double seal into and around the nose.

- **Easy Frame**

This low-profile frame is stable, durable and small, and ensures a clear line of sight. The delivery tube is lightweight, very flexible and designed to minimize pressure or pull on the face.

- **Headgear Options: Adjustable or Stretchwise™ Headgear**

a. **Adjustable Headgear** – This breathable, ergonomically designed headgear fits a greater range of head sizes and enables sizing adjustment to the patient's preference.







b. **Stretchwise Headgear** – This minimal headgear automatically adjusts to the wearer for a quick and easy fit.

Product Code	Description	Quantity	HCPCS Code
400421 F&P Pilairo Q	F&P Pilairo Q Nasal Pillows Mask (Complete System) includes: <ul style="list-style-type: none"> • Adjustable Headgear • Stretchwise Headgear • AirPillow Seal • Flexible Tube 	1 1 1 1	A7034 + A7035

Product Code	Description	HCPCS Code	
400HC125*	AirPillow Seal	A7033	
400HC328	Adjustable & Stretchwise Headgear (Combo Pack)**	A7035	
400HC326	Adjustable Headgear***	A7035	
400HC322	Stretchwise Headgear***	A9999	
400HC229	Elbow Cover & Diffuser		
400HC588	Headgear Clips & Buckle		
400HC224	Female Swivel		
400HC565*	Pilairo Q Nasal Pillows Mask (no headgear)		





† F&P Pilairo Q with Stretchwise Headgear
 * Reimbursable replacement parts
 ** Available in USA
 *** Not available in USA

Opus™ 360

Product Code	Description	HCPCS Code	
HC482	Opus 360 Complete system includes: <ul style="list-style-type: none"> • Mask Base and Flexitube Assembly • Headgear • Silicone Seal (Small, Medium and Large) • Male Swivel • Spare Tube-Anchoring Strap 	A7034 + A7035	
400HC116* 400HC117* 400HC118*	Silicone Pillow – Small Silicone Pillow – Medium Silicone Pillow – Large	A7033	
400HC310*	Headgear with Tube-Anchoring Strap (one size)	A7035	
400HC311	Tube-Anchoring Strap		
400HC312	Headgear – Tri-Glide Buckle		
400HC526*	Opus 360 Nasal Pillows Mask (no headgear)	A7034	

* Reimbursable replacement parts

Oracle™

Product Code	Description	HCPCS Code	
HC452	<p>Oracle</p> <p>Complete system includes:</p> <ul style="list-style-type: none"> • Flexitube Assembly • Headgear with Easy-Release Cord • Silicone Seal (Small and Large) • Storage Case 	A7044 + A9999	
400HC105 400HC106	<p>Silicone Seal – Small</p> <p>Silicone Seal – Large</p>		
400HC300	Headgear Accessory plus Clip	A9999	
900HC439	Exhaust Flow Diffuser Material – 10 Pack		
900HC440	Exhaust Flow Diffuser Cap plus Spare		
900HC441	Elbow Kit		
400HC204	Nasal Plugs – Small, Medium and Large		

Manufacturer 

Fisher & Paykel Healthcare Ltd
15 Maurice Paykel Place
East Tamaki, Auckland 2013

PO Box 14 348 Panmure
Auckland 1741
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Australia

Fisher & Paykel Healthcare
Pty Limited
36-40 New Street,
PO Box 167 Ringwood,
Melbourne, Victoria 3134,
Australia

Tel: +61 3 9879 5022
Fax: +61 3 9879 5232

Austria

Tel: 0800 29 31 23
Fax: 0800 29 31 22

Benelux

Tel: +31 40 216 3555
Fax: +31 40 216 3554

China

Tel: +86 20 3205 3486
Fax: +86 20 3205 2132

France

Tel: +33 1 6446 5201
Fax: +33 1 6446 5221

Germany

Tel: +49 7181 98599 0
Fax: +49 7181 98599 66

India

Tel: +91 80 4284 4000
Fax: +91 80 4123 6044

Irish Republic

Tel: 1800 409 011

Italy

Tel: +39 06 7839 2939
Fax: +39 06 7814 7709

Japan

Tel: +81 3 3661 7205
Fax: +81 3 3661 7206

Northern Ireland

Tel: 0800 132 189

Spain

Tel: +34 902 013 346
Fax: +34 902 013 379

Sweden

Tel: +46 8 564 76 680
Fax: +46 8 36 63 10

Switzerland

Tel: 0800 83 47 63
Fax: 0800 83 47 54

Taiwan

Tel: +886 2 8751 1739
Fax: +886 2 8751 5625

Turkey

Fisher Paykel Sağlık Ürünleri
Ticaret Limited Şirketi,
P.O. Box 06371 Ostim,
Alinteri Bulvarı 1161/1
Sokak
No. 12-14, Ankara, Turkey

Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK 

Fisher & Paykel Healthcare Ltd
Unit 16, Cordwallis Park
Clivemont Road,
Maidenhead, Berkshire
SL6 7BU, UK

Tel: +44 1628 626 136
Fax: +44 1628 626 146

USA/Canada

Tel: +1 800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001