









Product specifications

Oxygen concentrations ¹	At least 87% at all settings (maximum of 96%) Based on atmospheric pressure of 14.7 psia (101 kPa) at 70° F (21° C)	
Outlet pressure	20 psig	
Flow settings and pulse volumes	Setting Nom. pulse vol. at 20 BPM Max. minute vol. output 1 11 ml 220 ml/min 2 22 ml 440 ml/min 3 33 ml 660 ml/min 4 44 ml 880 ml/min 5 50 ml 1,000 ml/min +/- 25% or 6 ml, whichever is greater (average of 20 consecutive pulses) over the rated environmental range	
Sound level	43 dBA typical at settings 2 and 20 BPM; 46 dBA typical at settings 5 and 20 BPM	
Physical characteristics		
Size	With standard battery: 9.4 in. x 8.3 in. x 3.6 in. (23.9 cm x 21.1 cm x 9.1 cm) With extended battery: 10.2 in x 8.3 in x 2.6 in (25.9 cm x 21.1 cm x 9.1 cm)	
Weight	5 lbs. (2.3 kg) with standard battery installed 6 lbs. (2.7 kg) with extended battery installed	
User interface	Touch screen LCD	
Operating conditions		
Operating temperature	41° F to 95° F (5° C to 35° C)	
Relative humidity	15% to 93%	
Altitude	Up to 10,000 ft. (3,048 m)	
AC power		
Input	100 to 240 VAC, 50/60 Hz	
Output	19 VDC, 6.3 A max	
Power consumption	120 W while charging	
DC power		
Input	12 – 16 VDC	
Output	19 +/- 5% VDC, 6.3 A	
Audible alarm		
Sound level	56 dBA or greater when measured at 1 m	
Alarm type	No breath, low oxygen concentration, high breath rate, technical fault, low battery, no flow, depleted battery, external power failure	
Battery		
Chemistry	Lithium ion 14.4 VDC (nominal)	
Battery options	Standard battery, Extended battery (equals two standard batteries)	
Dimensions	Standard battery: 7.9 in. x 3.6 in. x 1.0 in. (20.1 cm x 9.1 cm x 2.5 cm) Extended battery: 7.9 in. x 3.6 in. x 1.9 in. (20.1 cm x 9.1 cm x 4.8 cm)	
Weight	Standard battery: 1.1 lbs. (0.50 kg) Extended battery 2.1 lbs. (0.95 kg)	
Watt-hour rating per standard battery	3.4 Ah x 3.6V x 8 cells = 97.9 Wh per standard battery	
Charge time	Standard battery: 4 hours maximum Extended battery: 8 hours maximum	
Typical battery duration on pulse s	etting 2, 20 BPM ²	
Standard battery	Up to 4.5 hours	
	Up to 9 hours fully charged batteries used with a new SimplyGo Mini System. Battery duration times will degrade conditions, use over time, and operational condition of the concentrator.	

Ordering information

Device	Part number
SimplyGo Mini system with standard battery includes: SimplyGo Mini portable oxygen concentrator, carrying case, one standard rechargeable lithium ion battery, DC and AC power cords, accessory bag, and user manual	US: 1113601 International: 1113604
SimplyGo Mini system with extended battery includes: SimplyGo Mini portable oxygen concentrator, carrying case, one extended rechargeable lithium ion battery, DC and AC power cords, accessory bag, and user manual	US: 1113602 International: 1113605



Replacement accessories	Part number
Standard lithium ion battery	1116816
Extended lithium ion battery	1116817
SimplyGo Mini AC power supply	1116818
SimplyGo Mini AC power cord	1082664
SimplyGo Mini DC power supply	1116819
SimplyGo Mini external battery charger	1116830
User Manual	1118724
Accessories - Bags - Brown	Part number
SimplyGo Mini carry bag and strap,	1119897

SimplyGo Mini carry bag and strap, brown	1119897
SimplyGo Mini accessory bag, brown	1116824
SimplyGo Mini backpack, brown	1116836
A DI I	Daut accordance

Accessories - Bags - Black	Part number
SimplyGo Mini carry bag and strap, black	1119928
SimplyGo Mini accessory bag, black	1119893
SimplyGo Mini backpack, black	1119894

Warranties	Part number
SimplyGo Mini 3-year system, 1-year sieve warranty	1120121
SimplyGo Mini 3-year system, 3-year sieve warranty	1122070
SimplyGo Mini 2-year sieve warranty	1120123

extension³

3-year sieve warranty extension is a post-sale warranty

³2-year sieve warranty extension is a post-sale warranty purchase option for those units purchased initially with the 3-year system, 1-year sieve warranty.

Learn more

For more information about SimplyGo Mini, call customer service at 1-800-345-6443.

Visit www.philips.us/simplygomini

The Right Fi

SimplyGo Mini is part of our Right Fit, a complete portfolio of oxygen products and programs inspired by patients and built for business.



© 2016 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Caution: US federal law restricts these devices to sale by or on the order of a physician.