

Made in France



HORIZON[®]
P5

Portable Oxygen Concentrator

Horizon[®] proposes premium and highly reliable range of oxygen concentrators designed for oxygen therapy patients who want to maintain their lifestyle.

Bluetooth connected

Light and intuitive

High O₂ purity

SCALEO
MEDICAL

BENEFITS

LIGHTWEIGHT
& COMPACT



20 SETTINGS FOR
PRECISE DOSAGE

from 0,25 to 5

EASY
MAINTENANCE



HORIZON[®]
P5



LONG LIFE
BATTERY

BLUETOOTH
CONNECTED*



3 YEARS DEVICE'S
WARRANTY**

* IF ACTIVATED

** BATTERY, CHARGER AND SIEVE BEDS WARRANTED FOR 1 YEAR

HORIZON[®]
P5 pack includes:

- Horizon[®] P5 portable concentrator
- Single battery or double battery
- Carry bag
- AC/DC power supply 100-240 V, 50-60 Hz
- User manual *not shown in picture*



Intensive use



Designed & manufactured in France

MOBILITY & ACCESSIBILITY

BACKPACK AS AN OPTION ADAPTED TO BE USED WITH THE CARRY BAG



Simple & intuitive



Easy to access filters



High oxygen purity



Charging and maintenance access

TECHNICAL FEATURES

Dimensions w x d x h	216 x 88 x 220 mm* 216 x 88 x 240 mm**	Power consumption	85 W
Weight	2,54 kg* and 2,90 kg**	Power compatibility	100 - 240 VAC / 50 HZ 60 HZ
Oxygen purity	> 87 %	Environmental tolerance	Temperature +5 to +40°C Relative humidity: 95 %
Oxygen flow rate setting	20 flow settings (from 0,25 to 5)	Storage environmental tolerance	Temperature: -25 to +70°C Relative humidity: 0 % - 95 %
Battery life	Up to 4,25 hours on setting 1* Up to 8,5 hours on setting 1**	Pulse flow in mL/min	Up to 1050 mL/min
Recharge time	Up to 2,5* hours and 4,5** hours	Sound level	37,3 dB at setting 2***
	*with simple battery		***as per 14-1 10/2018 MDS-HI

SCALEO
M E D I C A L



Designed and manufactured by
SCALEO Medical - Parc 2000, 107 rue Joe Dassin - 34 080 Montpellier - France
Tel. +33 467 72 94 86
www.scaleomedical.com - info@scaleomedical.com